

pNPP ELISA Substrate System (ELPN-500)

For Rapid Accurate Alkaline Phosphatase and ELISA Assays

DESCRIPTION

Para-nitrophenyl phosphate (pNPP) is a chromogenic substrate for most phosphatases such as alkaline phosphatases and acid phosphatases. The reaction yields *para*-nitrophenol, an intense yellow product that can be conveniently measured at 405 nm on a spectrophotometer. The pNPP ELISA Substrate System has been optimized for enzyme-linked immunosorbent assays (ELISA) which use alkaline phosphatase conjugated antibody.

APPLICATIONS

ELISA Assays with alkaline phosphatase conjugated antibody.

Characterization of Alkaline Phosphatase: determination of phosphatase activity and phosphatase mechanism of action.

Drug Discovery: high-throughput screen for inhibitors of alkaline phosphatase.

KIT CONTENTS

Catalog #	No. of Tests	Reagent	Buffer
ELPN-500	500	solid	100 mL

Storage conditions. Store the Reagent and Buffer at -20°C and at 4°C, respectively. Shelf life: 12 months.

This protocol can be downloaded online at www.bioassaysys.com.

Precautions: reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

PROCEDURES (96-WELL PLATE)

1. Reconstitute the Substrate system. Fresh reconstitution is recommended although the reconstituted Substrate may be stable for up to 4 weeks when stored at -20°C. Simply add 100 mL Buffer to the amber Reagent bottle. Mix thoroughly by inversion multiple times. After this is done, mark the bottle label as Reconstituted Substrate. Equilibrate to room temperature before use, for all subsequent procedures are carried out at room temperature.

2. After the secondary antibody (alkaline phosphatase conjugate) has been incubated and thoroughly washed to remove unbound conjugate, add 200 µL of the reconstituted pNPP substrate per well. Incubate the plate in the dark for approximately 30 min.

3. Read the plate at 405 nm on a plate reader.

If the plate can not be read immediately, add 50 µL Stop Solution (3 N NaOH, not provided) to terminate further color reaction. Mix by quickly tapping the plate.

GENERAL CONSIDERATIONS

Background too high. A blocking step is preferred. This can be done by using a 1-5% serum prior to the application of the primary antibody. Adequate washes (including 0.1% TWEEN) are necessary to remove unbound antibody.

Color too faint. Weak color development occurs when primary/or secondary antibody concentrations are too low. A control experiment can be run to determine the activity of the alkaline phosphatase conjugate. Alternatively, allow the reaction to proceed for a longer period of time or at a higher temperature (e.g. 37°C).

LITERATURE

Derango R, Page J (1996). The quantitation of coupled bead antibody by enzyme-linked immunosorbent assay. *J Immunoassay*. 17(2):145-53.

Brown-Augsburger P, Yue XM, Lockridge JA, McSwiggen JA, Kamboj D, Hillgren KM (2004). Development and validation of a sensitive, specific, and rapid hybridization-ELISA assay for determination of concentrations of a ribozyme in biological matrices. *J Pharm Biomed Anal*. 34(1):129-39.

Sirois M, Chartier C, Brun G, Hamel Y, Petitclerc C, Delage JM (1984). An enzyme-linked immunosorbent assay for the detection of complement components on red blood cells. *Am J Clin Pathol*. 82(1):67-73.

Huang Z, Olson NA, You W, Haugland RP (1992). A sensitive competitive ELISA for 2,4-dinitrophenol using 3,6-fluorescein diphosphate as a fluorogenic substrate. *J Immunol Methods*. 149:261-6.

Di Mauro S, Manes T, Hessle L, Kozlenkov A, Pizauro JM, Hoylaerts MF, Millan JL (2002). Kinetic characterization of hypophosphatasia mutations with physiological substrates. *J Bone Miner Res*. 17(8):1383-91.

Martinez MB, Flickinger M, Higgins L, Krick T, Nelsestuen GL (2001). Reduced outer membrane permeability of *Escherichia coli* O157:H7: suggested role of modified outer membrane porins and theoretical function in resistance to antimicrobial agents. *Biochemistry* 40(40):11965-74.